

Product datasheet for **MP212052**

Ptpn11 Mouse qPCR Primer Pair (NM_011202)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 19247 |
| Forward Sequence: | TGACTTCTGGCGGATGGTGTTTC |
| Reverse Sequence: | GACGTTCTAACACGCATGACC |
| ACCN: | BC057398 , BC059278 , NM_011202 , NM_011202.1 , NM_011202.2 , NM_011202.3 , BC003980 , BC028513 , BM899214 |
| UniProt ID: | P35235 |
| Synonyms: | 2700084A17Rik; AW536184; PTP1D; PTP2C; SAP-2; SH-PTP2; SH-PTP3; SHP-2; Shp2; Syp |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)