

Product datasheet for **MP212066**

Ptpn7 Mouse qPCR Primer Pair (NM_177081)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	320139
Forward Sequence:	GACCGCTACAAGACCATCTTGC
Reverse Sequence:	GGTGGCAATGTAGACCTTCTCC
ACCN:	<u>BC098475</u> , <u>NM_177081</u> , <u>NM_177081.1</u> , <u>NM_177081.2</u> , <u>NM_177081.3</u> , <u>BC138194</u> , <u>BC145190</u> , <u>NM_177081.4</u>
UniProt ID:	<u>Q8BUM3</u>
Synonyms:	BPTP-4; C920001D21Rik; HEPTP; LC-PTP; LPTP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)