

Product datasheet for **MP212087**

Ptprs Mouse qPCR Primer Pair (NM_011218)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19280
Forward Sequence:	TACCAGGTCCACTATGTGCGCA
Reverse Sequence:	AGTCTCAGGCTGGAGTTCGTT
ACCN:	<u>BC052462</u> , <u>BC083188</u> , <u>NM_011218</u> , <u>NM_011218.1</u> , <u>NM_011218.2</u> , <u>BC052462.1</u> , <u>BC054077</u>
UniProt ID:	<u>B0V2N1</u>
Synonyms:	AL022616; PTP; PTP-NU3; PTPNU-3; PTPsigma; Ptpt9; R-PTP-S; RPTPsigma
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)