

## Product datasheet for **MP212116**

### **Pcdh17 Mouse qPCR Primer Pair (NM\_001013753)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	219228
<b>Forward Sequence:</b>	CAAAGGCTCCTGCTGTGACATG
<b>Reverse Sequence:</b>	GACCAAGCACTCGGCATTCATC
<b>ACCN:</b>	<u><a href="#">BC089004</a>, <a href="#">NM_001013753</a>, <a href="#">NM_001013753.1</a>, <a href="#">NM_001013753.2</a>, <a href="#">BC079850</a></u>
<b>Synonyms:</b>	C030033F14Rik; Gm78
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
<b>Storage:</b>	Store at -20 °C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)