

## Product datasheet for **MP212225**

### **P2rx7 Mouse qPCR Primer Pair (NM\_011027)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	18439
<b>Forward Sequence:</b>	GAACACGGATGAGTCCTTCGTC
<b>Reverse Sequence:</b>	CAGTGCCGAAAACCAGGATGTC
<b>ACCN:</b>	<u><a href="#">NM_011027</a></u> , <u><a href="#">NM_011027.1</a></u> , <u><a href="#">NM_011027.2</a></u> , <u><a href="#">NM_011027.3</a></u> , <u><a href="#">BC062808</a></u> , <u><a href="#">BC128303</a></u> , <u><a href="#">BC128304</a></u> , <u><a href="#">BC141120</a></u> , <u><a href="#">BQ554327</a></u> , <u><a href="#">BY752747</a></u> , <u><a href="#">NM_011027.4</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q9Z1M0</a></u>
<b>Synonyms:</b>	AI467586; P2X(7); P2X7R
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)