

Product datasheet for **MP212274**

Qk Mouse qPCR Primer Pair (NM_021881)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19317
Forward Sequence:	AATGGTCCGAGGCAAAGGCTCA
Reverse Sequence:	GCTCTGTTCTGAGCATCTTCCAC
ACCN:	<u>BC053426</u> , <u>BC056346</u> , <u>NM_021881</u> , <u>NM_021881.1</u> , <u>NM_021881.2</u> , <u>BC053426.1</u> , <u>BB674938</u>
UniProt ID:	<u>Q9QYS9</u>
Synonyms:	I110003F05Rik; I500005P18; I(17)-1Wis; I17WisI; QkI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

