

Product datasheet for **MP212308**

Rab14 Mouse qPCR Primer Pair (NM_026697)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68365
Forward Sequence:	GCAGATTTGGGATACAGCAGGG
Reverse Sequence:	GGTGAGATTCCTTGCGTCTGTC
ACCN:	<u>BC009085</u> , <u>BC025139</u> , <u>BC056648</u> , <u>NM_026697</u> , <u>NM_026697.1</u> , <u>NM_026697.2</u> , <u>NM_026697.3</u> , <u>BU702057</u> , <u>NM_026697.4</u>
UniProt ID:	<u>Q91V41</u>
Synonyms:	0610030G24Rik; 2810475J17Rik; A830021G03Rik; AI314285; AI649155; D030017L14Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)