

Product datasheet for **MP212321**

Rab27a Mouse qPCR Primer Pair (NM_023635)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11891
Forward Sequence:	GAGCAAAGTTTCCTCAATGTCCG
Reverse Sequence:	CTTTCAGTCCCTCTGGTCTTC
ACCN:	<u>BC008173</u> , <u>BC009656</u> , <u>NM_023635</u> , <u>NM_023635.1</u> , <u>NM_023635.2</u> , <u>NM_023635.3</u> , <u>NM_023635.4</u> , <u>NM_023635.5</u> , <u>NM_023635.6</u> , <u>BB857060</u> , <u>BE948607</u> , <u>BP769437</u> , <u>BY151058</u>
UniProt ID:	<u>Q9ERI2</u>
Synonyms:	2210402C08Rik; 2410003M20Rik; 4933437C11Rik; ash
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)