

Product datasheet for **MP212331**

Rab34 Mouse qPCR Primer Pair (NM_033475)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19376
Forward Sequence:	ACAAGGCTACCATCGGAGTGGA
Reverse Sequence:	CGGTAGTAGGTGGAAGCAATGC
ACCN:	BC038638 , NM_033475 , NM_033475.1 , NM_033475.2 , NM_033475.3 , BB643883 , BC028481 , BY195955 , BY234941 , BY317641
UniProt ID:	Q64008
Synonyms:	AI326479; Narr; Rah; Rah1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)