

Product datasheet for **MP212349**

Rab5a Mouse qPCR Primer Pair (NM_025887)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	271457
Forward Sequence:	AAGCACAGTCCTATGCAGATGAC
Reverse Sequence:	GCTGAGTTGCACCAGGATTCTG
ACCN:	BC004842 , BC034370 , BC096481 , NM_025887 , NM_025887.1 , NM_025887.2 , NM_025887.3 , NM_025887.4 , BC004842.1
UniProt ID:	Q9CQD1
Synonyms:	2410015H04Rik; AI663973; AU021172; nnyRab5a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)