

Product datasheet for **MP212353**

Rab6a Mouse qPCR Primer Pair (NM_024287)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19346
Forward Sequence:	GGAAAGACGTCCTTGATCACCC
Reverse Sequence:	CTCTTGACCTGCTGTGCCAT
ACCN:	BC019118 , NM_024287 , NM_024287.1 , NM_024287.2 , NM_024287.3 , NM_024287.4 , BC026915 , BE134038 , BP765476
UniProt ID:	P35279
Synonyms:	2610028L11Rik; AA419671; Rab6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)