

Product datasheet for **MP212362**

Rabep1 Mouse qPCR Primer Pair (NM_019400)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54189
Forward Sequence:	CTGGAGGACTTCCTCAAGCAGA
Reverse Sequence:	ATCTGCTCCTGAACGTCCCTCT
ACCN:	<u>BC129853</u> , <u>BC129854</u> , <u>NM_019400</u> , <u>NM_019400.1</u> , <u>NM_019400.2</u> , <u>NM_019400.3</u> , <u>BC003921</u> , <u>BC009814</u> , <u>BC060166</u> , <u>BC080274</u> , <u>BY133091</u> , <u>BY718394</u>
UniProt ID:	<u>O35551</u>
Synonyms:	Rab5ep; rabaptin-5; rabaptin-5alpha
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)