

Product datasheet for **MP212366**

Racgap1 Mouse qPCR Primer Pair (NM_012025)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26934
Forward Sequence:	GGAGCCAAGAAGTACACAGAC
Reverse Sequence:	ACTTGATCCGCTTCCACACGG
ACCN:	<u>BC010715</u> , <u>NM_012025</u> , <u>NM_012025.1</u> , <u>NM_012025.2</u> , <u>NM_012025.3</u> , <u>NM_012025.4</u> , <u>NM_012025.5</u> , <u>NM_012025.6</u> , <u>NM_012025.7</u> , <u>BY309801</u> , <u>BY439051</u>
UniProt ID:	<u>Q9WVM1</u>
Synonyms:	AI227039; AI327394; Band25; gtl11; GTPase; MgcRacGAP; mKIAA1478
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)