

## Product datasheet for **MP212375**

### Rad23b Mouse qPCR Primer Pair (NM\_009011)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19359
Forward Sequence:	GCCTTGCTACAGCAGATAGGTC
Reverse Sequence:	CTCCACCTTGACTGCCTGCTTC
ACCN:	<u><a href="#">BC027747</a></u> , <u><a href="#">BC068193</a></u> , <u><a href="#">NM_009011</a></u> , <u><a href="#">NM_009011.1</a></u> , <u><a href="#">NM_009011.2</a></u> , <u><a href="#">NM_009011.3</a></u> , <u><a href="#">NM_009011.4</a></u> , <u><a href="#">BE989853</a></u> , <u><a href="#">BY140385</a></u>
UniProt ID:	<u><a href="#">P54728</a></u>
Synonyms:	0610007D13Rik; AV001138; HR23B; mHR23B; p58
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)