

Product datasheet for **MP212386**

Rad54l2 Mouse qPCR Primer Pair (NM_030730)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81000
Forward Sequence:	TGATCCTGGTGCGGCTTTCTCA
Reverse Sequence:	CAGCACATCAGGATGGTTCCAG
ACCN:	BC133714 , NM_030730 , NM_030730.1 , NM_030730.2 , BF455201
UniProt ID:	Q99NG0
Synonyms:	Arip4; D130058C05; G630026H09Rik; Srisnf2l
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)