

## Product datasheet for **MP212432**

### Rora Mouse qPCR Primer Pair (NM\_013646)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19883
Forward Sequence:	CAGAGCAATGCCACCTACTCCT
Reverse Sequence:	CTGCTTCTGGACATCCGACCA
ACCN:	<u><a href="#">BC003757</a></u> , <u><a href="#">NM_013646</a></u> , <u><a href="#">NM_013646.1</a></u> , <u><a href="#">NM_013646.2</a></u>
UniProt ID:	<u><a href="#">P51448</a></u>
Synonyms:	9530021D13Rik; nmf267; Nr1f1; ROR1; ROR2; ROR3; sg; staggerer; tmgc26
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)