

Product datasheet for **MP212460**

Rhoj Mouse qPCR Primer Pair (NM_023275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80837
Forward Sequence:	GCCTTATGTGCTCATCGGAACC
Reverse Sequence:	CTTCACGCCGTGCTCATAGGTG
ACCN:	BC043719 , BC060193 , NM_023275 , NM_023275.1 , NM_023275.2 , BC024626 , BC054464 , BY004626
UniProt ID:	Q9ER71
Synonyms:	1110005O19Rik; Arhj; AW210585; TC10L; TCL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

