

Product datasheet for **MP212530**

Ror1 Mouse qPCR Primer Pair (NM_013845)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26563
Forward Sequence:	GACCTGTGTGACATCCCAGCAT
Reverse Sequence:	AGAAGAGGAAGGCGATAGCCAG
ACCN:	<u>NM_013845</u> , <u>NM_013845.1</u> , <u>NM_013845.2</u> , <u>NM_013845.3</u> , <u>NM_013845.4</u> , <u>NM_013845.5</u> , <u>BC138794</u> , <u>BB652069</u> , <u>BC095933</u> , <u>BC138810</u>
UniProt ID:	<u>Q9Z139</u>
Synonyms:	2810404D04Rik; Ntr; Ntrkr1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)