

Product datasheet for **MP212533**

Ripk3 Mouse qPCR Primer Pair (NM_019955)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56532
Forward Sequence:	GAAGACACGGCACTCCTTGTA
Reverse Sequence:	CTTGAGGCAGTAGTTCTTGGTGG
ACCN:	BC029210 , NM_019955 , NM_019955.1 , NM_019955.2 , BF658559 , BY206463 , BY209344
UniProt ID:	Q9QZL0
Synonyms:	2610528K09Rik; AW107945; Rip3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)