

## Product datasheet for **MP212539**

### **Rbpj Mouse qPCR Primer Pair (NM\_009035)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	19664
<b>Forward Sequence:</b>	TGGCTACATCCATTACGGGCAG
<b>Reverse Sequence:</b>	GTGGAGTTGTGATACAGGGTCG
<b>ACCN:</b>	<u><a href="#">BC051387</a></u> , <u><a href="#">NM_009035</a></u> , <u><a href="#">NM_009035.1</a></u> , <u><a href="#">NM_009035.3</a></u> , <u><a href="#">NM_009035.4</a></u> , <u><a href="#">NM_009035.5</a></u> , <u><a href="#">BB839243</a></u> , <u><a href="#">BC035299</a></u> , <u><a href="#">BY107008</a></u> , <u><a href="#">BY155397</a></u> , <u><a href="#">BY171919</a></u> , <u><a href="#">BY716845</a></u>
<b>UniProt ID:</b>	<u><a href="#">P31266</a></u>
<b>Synonyms:</b>	AI843960; CBF1; Igkjrb; Igkrsbp; RBP-J; RBP-J kappa; RBP-Jkappa; RBPjk; Rbpsuh
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)