

Product datasheet for **MP212604**

Rxfp4 Mouse qPCR Primer Pair (NM_181817)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	242093
Forward Sequence:	AGGTAACGTGGTCAGCGTGTG
Reverse Sequence:	CCAGTGGAAGTCCATTGCTGAC
ACCN:	<u>BC107269</u> , <u>NM_181817</u> , <u>NM_181817.1</u> , <u>BC107270</u>
UniProt ID:	<u>Q7TQP4</u>
Synonyms:	GPCR142; Gpr100; Rln3r2; SALPR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)