

## Product datasheet for **MP212718**

### **Rbp3 Mouse qPCR Primer Pair (NM\_015745)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	19661
<b>Forward Sequence:</b>	CACCTCCCATTCTGCTGGACAA
<b>Reverse Sequence:</b>	ACACCACTGGTCAGGATGGCTA
<b>ACCN:</b>	<a href="#">BC017610</a> , <a href="#">BC042480</a> , <a href="#">BC083133</a> , <a href="#">NM_015745</a> , <a href="#">NM_015745.1</a> , <a href="#">NM_015745.2</a> , <a href="#">BC023069</a> , <a href="#">BC034121</a> , <a href="#">BC035273</a> , <a href="#">BC040282</a> , <a href="#">BC042509</a> , <a href="#">NM_015745.3</a>
<b>UniProt ID:</b>	<a href="#">P49194</a>
<b>Synonyms:</b>	lrbp; Rbp-3
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)