

Product datasheet for **MP212722**

Rdh10 Mouse qPCR Primer Pair (NM_133832)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	98711
Forward Sequence:	TGGTCAACTGCCACGCACACTT
Reverse Sequence:	CCTCAACTCCAGCAGTGCTGAA
ACCN:	<u>BC019796</u> , <u>NM_133832</u> , <u>NM_133832.1</u> , <u>NM_133832.2</u> , <u>NM_133832.3</u> , <u>BC038334</u> , <u>BC052546</u>
UniProt ID:	<u>Q8VCH7</u>
Synonyms:	3110069K09Rik; 4921506A21Rik; AI875664; AW549993; D1ErtD762e; m366Asp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)