

Product datasheet for **MP212723**

Rdh11 Mouse qPCR Primer Pair (NM_021557)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17252
Forward Sequence:	AGTCACTGGTGCTAACACAGGC
Reverse Sequence:	TGTTCCCTGTGACGGCTTGGAT
ACCN:	<u>BC018261</u> , <u>NM_021557</u> , <u>NM_021557.1</u> , <u>NM_021557.2</u> , <u>NM_021557.3</u> , <u>NM_021557.4</u> , <u>NM_021557.5</u> , <u>BC036273</u> , <u>NM_021557.6</u>
UniProt ID:	<u>Q9QYF1</u>
Synonyms:	2610319N22Rik; AI428145; Arsd1; AU045252; C85936; CGI-82; HCBP12; M42C60; Mdt1; Psdr1; SCALD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)