

## Product datasheet for **MP212755**

### Rnd1 Mouse qPCR Primer Pair (NM\_172612)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                              |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| Gene ID:          | 223881                                                                                                                                                                                                                                                                                                                                       |
| Forward Sequence: | CCTATGTGCCTACCGTGTTGAG                                                                                                                                                                                                                                                                                                                       |
| Reverse Sequence: | CAGCAATACCGCATCTGAGTCG                                                                                                                                                                                                                                                                                                                       |
| ACCN:             | <u><a href="#">BC048531</a>, <a href="#">NM_172612</a>, <a href="#">NM_172612.1</a>, <a href="#">NM_172612.2</a>, <a href="#">NM_172612.3</a>, <a href="#">BC048531.1</a></u>                                                                                                                                                                |
| UniProt ID:       | <u><a href="#">Q8BLR7</a></u>                                                                                                                                                                                                                                                                                                                |
| Synonyms:         | A830014L09Rik; Arhs                                                                                                                                                                                                                                                                                                                          |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                              |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                 |



[View online »](#)