

Product datasheet for **MP212758**

Arhgdia Mouse qPCR Primer Pair (NM_133796)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	192662
Forward Sequence:	GACAAGGACGATGAAAGCCTCC
Reverse Sequence:	CCTGTCAGGTCCAGTTCCAGAG
ACCN:	<u>BC004732, BC086755, NM_133796, NM_133796.2, NM_133796.3, NM_133796.4, NM_133796.5, NM_133796.6, NM_133796.7, NM_133796.8</u>
UniProt ID:	<u>Q99PT1</u>
Synonyms:	5330430M07Rik; Gdi-1; Rho-GDI; RhoDGI; RhoGDI-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)