

Product datasheet for **MP212793**

Arhgef4 Mouse qPCR Primer Pair (NM_183019)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226970
Forward Sequence:	TCAAGCCAGGAGTCAGCAGAGA
Reverse Sequence:	CATCTACCACCTCCAGGTCATC
ACCN:	<u>NM_183019</u> , <u>NM_183019.1</u> , <u>NM_183019.2</u> , <u>BC139397</u> , <u>BB609907</u> , <u>BC055804</u> , <u>BC094608</u> , <u>BC139371</u> , <u>NM_183019.3</u>
UniProt ID:	<u>Q7TNR9</u>
Synonyms:	9330140K16Rik; Asef; ENSMUSG00000070955
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)