

Product datasheet for **MP212797**

Arhgef2 Mouse qPCR Primer Pair (NM_008487)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16800
Forward Sequence:	GATGTGCTCGTGTCTCCAGG
Reverse Sequence:	ACATCTCAGGCGGTCCAGAACT
ACCN:	<u>BC006589</u> , <u>NM_008487</u> , <u>NM_008487.1</u> , <u>NM_008487.2</u> , <u>NM_008487.3</u> , <u>BX521217</u> , <u>BY316631</u>
UniProt ID:	<u>Q60875</u>
Synonyms:	AA408978; GEF; GEF-H1; GEFH1; Lbcl1; Lfc; LFP40; mKIAA0651; P40
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)