

## Product datasheet for **MP212801**

### Arhgef28 Mouse qPCR Primer Pair (NM\_012026)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	110596
Forward Sequence:	GCTCAGGAGTTTGAAGCAGAGTC
Reverse Sequence:	TGAGCAGAGTCTGGATGTGGTG
ACCN:	<u><a href="#">NM_012026</a></u> , <u><a href="#">NM_012026.1</a></u> , <u><a href="#">NM_012026.2</a></u> , <u><a href="#">BC138030</a></u> , <u><a href="#">BC022978</a></u> , <u><a href="#">BC057883</a></u> , <u><a href="#">BC075734</a></u> , <u><a href="#">BC098331</a></u> , <u><a href="#">BY778560</a></u>
UniProt ID:	<u><a href="#">P97433</a></u>
Synonyms:	9230110L08Rik; AI323540; D13Bwg1089e; p190RhoGEF; Rgnef; RhoGEF; Rhoip2; RIP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)