

Product datasheet for **MP212859**

Rpl14 Mouse qPCR Primer Pair (NM_025974)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67115
Forward Sequence:	CGCCAGAAGTATGTCCGAAAAGC
Reverse Sequence:	CGATCAAAGTCTGTCATCTTGGC
ACCN:	BC092249 , NM_025974 , NM_025974.1 , NM_025974.2 , BC013249 , BC013251 , BC037202 , BC042937
UniProt ID:	Q9CR57
Synonyms:	3100001N19Rik; AA407502; AL022816; CAG-ISL-7; CTG-B33; L14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)