

Product datasheet for **MP212861**

Rpl17 Mouse qPCR Primer Pair (NM_001002239)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319195
Forward Sequence:	CGCGCGTATAATGGTGGAGTTG
Reverse Sequence:	TAAGTTCAGCGTTACTCTCAGCAT
ACCN:	BC052940 , BC054424 , BC090990 , BC091759 , BC092091 , BC106165 , NM_001002239 , NM_001002239.1 , NM_001002239.2 , NM_001002239.3 , BC003896 , BC090401 , BC094640 , BY031981
Synonyms:	MGC103289; MGC118163
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)