

## Product datasheet for **MP212865**

### Rpl21 Mouse qPCR Primer Pair (NM\_019647)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19933
Forward Sequence:	GAGCAGAGACAGCTTCTCTGAAG
Reverse Sequence:	CCTCACAAGTGTGCTTCTCTGG
ACCN:	<u><a href="#">BC086904</a></u> , <u><a href="#">BC086905</a></u> , <u><a href="#">BC086935</a></u> , <u><a href="#">BC089582</a></u> , <u><a href="#">BC090256</a></u> , <u><a href="#">BC094410</a></u> , <u><a href="#">BC096404</a></u> , <u><a href="#">BC096572</a></u> , <u><a href="#">BC106116</a></u> , <u><a href="#">NM_019647</a></u> , <u><a href="#">NM_019647.1</a></u> , <u><a href="#">NM_019647.2</a></u> , <u><a href="#">NM_019647.3</a></u> , <u><a href="#">NM_019647.4</a></u> , <u><a href="#">NM_019647.5</a></u> , <u><a href="#">NM_019647.6</a></u> , <u><a href="#">BC096039</a></u> , <u><a href="#">BC138299</a></u> , <u><a href="#">BC138300</a></u> , <u><a href="#">BQ031254</a></u> , <u><a href="#">BY264537</a></u>
Synonyms:	8430440E03Rik; L21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)