

Product datasheet for **MP212876**

Rpl3 Mouse qPCR Primer Pair (NM_013762)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27367
Forward Sequence:	CAGCGATGAGTGTAAAAGGCGC
Reverse Sequence:	GCTGCTTCTTGCCTGTGCATC
ACCN:	BC009655 , BC083134 , BC094059 , NM_013762 , NM_013762.1 , NM_013762.2 , BC085495 , BC098497 , BC158039 , BY581903
UniProt ID:	P27659
Synonyms:	F2; J1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)