

Product datasheet for **MP212887**

Rpl37a Mouse qPCR Primer Pair (NM_009084)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19981
Forward Sequence:	AATCAGCCAGCACGCCAAGTAC
Reverse Sequence:	GCCACTGTTTTTCATGCAGGAACC
ACCN:	<u>BC086796</u> , <u>BC096413</u> , <u>NM_009084</u> , <u>NM_009084.1</u> , <u>NM_009084.2</u> , <u>NM_009084.3</u> , <u>NM_009084.4</u> , <u>BP773503</u>
UniProt ID:	<u>P61514</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)