

Product datasheet for **MP212896**

Rpl7 Mouse qPCR Primer Pair (NM_011291)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19989
Forward Sequence:	TCGCAGAGTTGAAGGTGAAGCG
Reverse Sequence:	CCATCCGAATCTCAGTGCGGTA
ACCN:	BC025909 , BC051261 , BC086786 , BC096452 , NM_011291 , NM_011291.1 , NM_011291.2 , NM_011291.3 , NM_011291.4 , NM_011291.5 , BB839528 , BY441493
UniProt ID:	P14148
Synonyms:	Rpl7a; Surf-3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)