

Product datasheet for **MP212900**

Rpl9 Mouse qPCR Primer Pair (NM_011292)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20005
Forward Sequence:	CATCAACGTGGAGCTGAGTCTTC
Reverse Sequence:	GACCGCATCTTGTATCGGAAGC
ACCN:	BC013165 , BC081435 , BC083166 , BC083329 , BC086937 , BC089319 , BC093526 , NM_011292 , NM_011292.1 , NM_011292.2 , BC069963 , BY072215
UniProt ID:	P51410
Synonyms:	MGC102393
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)