

Product datasheet for **MP212903**

Rps12 Mouse qPCR Primer Pair (NM_011295)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20042
Forward Sequence:	AAGACCGCCCTCATCCACGATG
Reverse Sequence:	TTGGTGCTCAGCAGAAAGTGCC
ACCN:	<u>BC002079</u> , <u>BC089338</u> , <u>BC089339</u> , <u>BC090257</u> , <u>BC092044</u> , <u>BC096653</u> , <u>BC099377</u> , <u>NM_011295</u> , <u>NM_011295.1</u> , <u>NM_011295.2</u> , <u>NM_011295.3</u> , <u>NM_011295.4</u> , <u>NM_011295.5</u> , <u>NM_011295.6</u> , <u>BY128264</u>
Synonyms:	MGC102111; MGC102112; MGC102499; MGC117504
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)