

Product datasheet for **MP212907**

Rps15a Mouse qPCR Primer Pair (NM_170669)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	267019
Forward Sequence:	ACCTCACAGGAAGGTTGAACAAG
Reverse Sequence:	TAATGCCAGCCGAGGTTGTCAG
ACCN:	BC054792 , BC055697 , BC076563 , BC081428 , BC086885 , BC087867 , NM_170669 , NM_170669.1 , NM_170669.2 , BC051205
UniProt ID:	P62245
Synonyms:	A630031B11Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)