

## Product datasheet for **MP212932**

### **Rps6ka3 Mouse qPCR Primer Pair (NM\_148945)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	110651
<b>Forward Sequence:</b>	TAACCGCAGAGGTCACACTCAG
<b>Reverse Sequence:</b>	CTCAGAAACTGTGGCATCCCGA
<b>ACCN:</b>	<a href="#">NM_148945</a> , <a href="#">NM_148945.1</a> , <a href="#">NM_148945.2</a> , <a href="#">BC008162</a> , <a href="#">BC038683</a> , <a href="#">BC150156</a> , <a href="#">BC150478</a> , <a href="#">BC165985</a>
<b>UniProt ID:</b>	<a href="#">P18654</a>
<b>Synonyms:</b>	MAPKAPK-1b; MPK-9; p90RSK3; pp90RSK2; Rsk2; S6K-alpha3
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)