

Product datasheet for **MP213041**

Tyw5 Mouse qPCR Primer Pair (NM_001037742)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68736
Forward Sequence:	CCTGAGTCAAGTTGGAGGAACG
Reverse Sequence:	GTTTCTTCGGCTGCTCTCTGGA
ACCN:	<u>BC132289</u> , <u>NM_001037742</u> , <u>NM_001037742.1</u> , <u>NM_001037742.2</u> , <u>BC023868</u> , <u>BC058590</u> , <u>BC099495</u> , <u>BC144765</u> , <u>BC144766</u> , <u>BU697704</u>
UniProt ID:	<u>A2RSX7</u>
Synonyms:	1110034B05Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)