

Product datasheet for **MP213056**

Lgalsl Mouse qPCR Primer Pair (NM_173752)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216551
Forward Sequence:	CCAGTTCAAGCCGACGTGTACT
Reverse Sequence:	TCCACAGGTCAAGCTGATGGCA
ACCN:	<u>BC019131</u> , <u>NM_173752</u> , <u>NM_173752.1</u> , <u>NM_173752.2</u> , <u>NM_173752.3</u> , <u>NM_173752.4</u> , <u>BC005536</u> , <u>BE979710</u>
UniProt ID:	<u>Q8VED9</u>
Synonyms:	1110067D22Rik; A530071M23; Grpa; Hspc159; Lgalsla
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)