

## Product datasheet for **MP213188**

### Mtfp1 Mouse qPCR Primer Pair (NM\_026443)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67900
Forward Sequence:	GCTATGGTGTGTCCAGCTCCTA
Reverse Sequence:	AGACTGCCACACAAAGGTGTCG
ACCN:	<a href="#">BC019205</a> , <a href="#">BC058737</a> , <a href="#">NM_026443</a> , <a href="#">NM_026443.1</a> , <a href="#">NM_026443.2</a> , <a href="#">NM_026443.3</a> , <a href="#">NM_026443.4</a> , <a href="#">BC012242</a>
UniProt ID:	<a href="#">Q9CRB8</a>
Synonyms:	1700020C11Rik; 2610507A21Rik; AV005788; Mtp18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)