

Product datasheet for **MP213206**

1700025G04Rik Mouse qPCR Primer Pair (NM_197990)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69399
Forward Sequence:	CCAGAATGGAGATGTGTTTGCCG
Reverse Sequence:	CCTGTGTTTCGAGCTGGCACTT
ACCN:	<u>BC022177</u> , <u>BC034723</u> , <u>BC055912</u> , <u>NM_197990</u> , <u>NM_197990.1</u> , <u>NM_197990.2</u> , <u>NM_197990.3</u> , <u>BB840268</u>
UniProt ID:	<u>Q8K207</u>
Synonyms:	2610510E02Rik; 2900042O04Rik; ENSMUSG00000073540
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)