

Product datasheet for **MP213599**

Mvb12b Mouse qPCR Primer Pair (NM_175184)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72543
Forward Sequence:	CCAGTCATCAGCCTCTTCTACC
Reverse Sequence:	GCTGGTACTCATAGTCTGTGCTG
ACCN:	<u>BC049129</u> , <u>BC059907</u> , <u>NM_175184</u> , <u>NM_175184.1</u> , <u>NM_175184.2</u> , <u>NM_175184.3</u> , <u>NM_175184.4</u> , <u>BC051542</u>
UniProt ID:	<u>Q6KAU4</u>
Synonyms:	2610200O14Rik; 2610528K11Rik; AI414895; Fam125b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)