

Product datasheet for **MP213642**

Cyb5rl Mouse qPCR Primer Pair (NM_175471)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	230582
Forward Sequence:	TGGAGAAGGTCACCAAGGACAC
Reverse Sequence:	CTGATGGGAGTGTAGGCTCTCT
ACCN:	<u>BC043687, NM_175471, NM_175471.1, NM_175471.2, NM_175471.3, BC004750, NM_175471.4</u>
UniProt ID:	<u>B1AS42</u>
Synonyms:	2810410C14Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)