

## Product datasheet for **MP213672**

### **Myl12a Mouse qPCR Primer Pair (NM\_026064)**

#### **Product data:**

|                          |  |
|--------------------------|--|
| <b>Product Type:</b>     | qPCR Primer Pairs  |
| <b>Gene ID:</b>          | 67268  |
| <b>Forward Sequence:</b> | GTTCCGCATGTTTGACCAGTCC   |
| <b>Reverse Sequence:</b> | CGTCCAGGTATTCGTCAGTTGG   |
| <b>ACCN:</b>             | <u><a href="#">BC010265</a></u> , <u><a href="#">BC029180</a></u> , <u><a href="#">NM_026064</a></u> , <u><a href="#">NM_026064.1</a></u> , <u><a href="#">NM_026064.2</a></u> , <u><a href="#">BC006719</a></u> , <u><a href="#">BC020282</a></u> , <u><a href="#">BC089050</a></u> , <u><a href="#">BY258797</a></u> , <u><a href="#">BY319676</a></u> |
| <b>Synonyms:</b>         | 2900073G15Rik; AI838292; C80200  |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.             |
| <b>Storage:</b>          | Store at -20°C.  |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)