

Product datasheet for **MP213680**

3110001I22Rik Mouse qPCR Primer Pair (NM_025653)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66598
Forward Sequence:	CCATTACTCAGAACAGCAGTCTC
Reverse Sequence:	GGAGAATGGCTTCCACTGGATC
ACCN:	BC038880 , NM_025653.1 , NM_025653.2
Synonyms:	3010001A07Rik; 3110001I22Rik; AI666707; AW107665; Bar; Rnf47
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)