

Product datasheet for **MP213715**

Msl1 Mouse qPCR Primer Pair (NM_028722)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74026
Forward Sequence:	GATAGAGGAGCTGAAGTCCGAG
Reverse Sequence:	TGGAGCAGCTTGTGCCTTTCCT
ACCN:	BC043039 , BC055715 , BC058629 , NM_028722 , NM_028722.1 , NM_028722.2 , BM946681 , BQ176491 , NM_028722.3
UniProt ID:	Q6PDM1
Synonyms:	2810017F12Rik; 4121402D02Rik; 4930463F05Rik; AA682082; Msl-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)