

Product datasheet for **MP213722**

Chst15 Mouse qPCR Primer Pair (NM_029935)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	77590
Forward Sequence:	GTGTGGAACCACTGACCTGTAC
Reverse Sequence:	ATCTTCCACAGGGTAGCGGTCT
ACCN:	BC031443 , NM_029935 , NM_029935.1 , NM_029935.2 , NM_029935.3 , NM_029935.4 , NM_029935.5 , NM_029935.6
UniProt ID:	Q91XQ5
Synonyms:	4631426J05Rik; BRAG; GalNAcS-6ST; MAd5; mKIAA0598
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)