

## Product datasheet for **MP213733**

### Vsir Mouse qPCR Primer Pair (NM\_028732)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 74048  |
| Forward Sequence: | AAACCACCACCCAGAACAACGG   |
| Reverse Sequence: | ACTGTCCTGCTCATTAGACGCC   |
| ACCN:             | <u><a href="#">BC003967</a></u> , <u><a href="#">NM_028732</a></u> , <u><a href="#">NM_028732.1</a></u> , <u><a href="#">NM_028732.2</a></u> , <u><a href="#">NM_028732.3</a></u> , <u><a href="#">NM_028732.4</a></u> , <u><a href="#">BY202363</a></u>   |
| UniProt ID:       | <u><a href="#">Q9D659</a></u>  |
| Synonyms:         | 4632428N05Rik; Dies1; PD-1H; VISTA   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)